

The  
United  
States  
of  
America



**The Director of the United States  
Patent and Trademark Office**

*Has received an application for a patent for a new and useful invention. The title and description of the invention are enclosed. The requirements of law have been complied with, and it has been determined that a patent on the invention shall be granted under the law.*

*Therefore, this*

**United States Patent**

*Grants to the person(s) having title to this patent the right to exclude others from making, using, offering for sale, or selling the invention throughout the United States of America or importing the invention into the United States of America for the term set forth below, subject to the payment of maintenance fees as provided by law.*

*If this application was filed prior to June 8, 1995, the term of this patent is the longer of seventeen years from the date of grant of this patent or twenty years from the earliest effective U.S. filing date of the application, subject to any statutory extension.*

*If this application was filed on or after June 8, 1995, the term of this patent is twenty years from the U.S. filing date, subject to any statutory extension. If the application contains a specific reference to an earlier filed application or applications under 35 U.S.C. 120, 121 or 365(c), the term of the patent is twenty years from the date on which the earliest application was filed, subject to any statutory extensions.*

*Jon W. I. Dudas*

*Director of the United States Patent and Trademark Office*



US007276040B2

(12) **United States Patent**  
**Gomez-de-Diego**

(10) **Patent No.:** **US 7,276,040 B2**  
(45) **Date of Patent:** **Oct. 2, 2007**

(54) **TRACTION APPARATUS FOR THE PENIS**

(76) **Inventor:** **Eduardo Gomez-de-Diego, C/**  
**Princessa 22, 6 cha, Madrid (ES) 28008**

(\*) **Notice:** Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 510 days.

(21) **Appl. No.:** **10/829,006**

(22) **Filed:** **Apr. 21, 2004**

(65) **Prior Publication Data**  
US 2004/0215055 A1 Oct. 28, 2004

(30) **Foreign Application Priority Data**  
Apr. 28, 2003 (EP) ..... 03380101

(51) **Int. Cl.**  
**A61F 5/00** (2006.01)

(52) **U.S. Cl.** ..... **602/36; 600/39**

(58) **Field of Classification Search** ..... **602/36;**  
**600/38, 39, 41; 128/842**  
See application file for complete search history.

(56) **References Cited**

**FOREIGN PATENT DOCUMENTS**

WO WO9626691 \* 9/1996

\* cited by examiner

*Primary Examiner*—Kim M. Lewis

(74) *Attorney, Agent, or Firm*—Allen, Dyer, Doppelt, Milbrath & Gilchrist, P.A.

(57) **ABSTRACT**

An apparatus and method for applying traction to a patient's penis includes at least one base member having a generally arcuate or circular periphery of sufficient dimension for the patient to insert the penis in the apparatus. At least one base member connector is positioned on the base member. A plurality of extension members is connected with the base member connector and extends distally therefrom so that at least one support member connected with the plurality of extension members is distally spaced apart from the base member, the support member having a surface for thereon supporting the penis. A retaining collar connected to the support member adjacent the support surface is positioned so as to at least partially encircle the shaft of the penis thereby securing the distal penis to the support member.

**37 Claims, 4 Drawing Sheets**

